


Basics of X-ray Scattering

Manfred Roessle
Luebeck University of Applied Science

X-Ray scattering


The Electromagnetic Spectrum



X-ray Scattering

Description of a plane electro-magnetic wave

Electromagnetic waves:




The electromagnetic waves that compose electromagnetic radiation can be imagined as a self-propagating transverse oscillating wave of electric and magnetic fields. This diagram shows a plane linearly polarized EMR wave propagating from left to right. The **electric field** is in a vertical plane and the **magnetic field** in a horizontal plane. The two types of fields are always in phase with each other, and no matter how powerful, have a ratio of electric to magnetic intensity which is fixed and never varies.

Source: Wikipedia


X-Ray scattering

X-ray scattering



X-rays are scattered at the electrons of the atomic shell. During the scattering process the electron starts oscillating. It becomes a dipole and a spherical wave is sent out. The wavelength and energy of the scattered wave does not change (elastic scattering).

X-ray scattering Geometrical approach



Scattering vector $\vec{s} = (\vec{k} - \vec{k}_0)$

Scattering angle 2Θ

X-ray scattering

Scattering from a single electron

For a single scattering process the amplitude A_j of scattered X-ray photons can be described as a plain wave scattered by an ensemble of atoms:

$$A_j = b_j e^{-i \frac{2\pi}{\lambda} (\vec{k}_0 - \vec{k}) \cdot \vec{r}}$$

With b_j is the scattering cross section, the \vec{r}_j describes the inner distance vector and the vector \vec{k}_0 is the wave vector of the incident wave.

X-ray scattering

Scattering from an ensemble

Within an particle, the scattering amplitudes of all atoms have to be summed up:

$$A(\vec{k}) = \sum_j b_j e^{-i \frac{2\pi}{\lambda} (\vec{k}_0 - \vec{k}) \cdot \vec{r}_j}$$

Or, using $\vec{s} = (\vec{k} - \vec{k}_0)$ as **scattering vector**

$$A(\vec{s}) = \sum_j b_j e^{-i \frac{2\pi}{\lambda} \vec{s} \cdot \vec{r}_j}$$

X-ray scattering

Scattering intensity

The scattering amplitude is experimentally not accessible, but the scattering intensity. The intensity is the product of the scattering amplitude with its complex conjugate and results to:

$$AA^* = I(\vec{s}) = \sum_j \sum_k b_j b_k e^{\vec{s} \cdot \vec{r}_{jk}}$$

With the vector \vec{r}_{jk} as inner distance vector, pointing from the j^{th} scattering center to the k^{th} scattering center.


The b_j are the atomic scattering densities. They depend on the number of electrons of the atom and the heavier an atom the higher is b_j .

X-ray scattering

Atomic scattering factors

The atomic scattering factor is the Fourier Transform of the atomic electron density.

It can be calculated by quantum mechanics and can be found at International Tables of Crystallography, Volume C



X-ray scattering

Electron density

As seen on the atomic scattering factors, the scattering density of the depends not only on the number of electrons, but as well on the distribution of these electrons around the atom.

This electron density ρ can be defined as follows:

$$b_j = \rho_j dV_j$$

With the dV_j as the volume element surrounding the scattering center.

The reconstitution of the electron density $\rho(r)$ and its distribution in the particle is the aim of every X-ray scattering or diffraction experiment.

X-Ray scattering

X-ray scattering

Electron density

The effective scattering density depends on the electron density of the surrounding of the atom.

A net electron density can be defined by subtracting the electron density of the environment ρ_S :

$$\Delta\rho(\vec{r}) = \rho(\vec{r}) - \rho_S$$

The average value of the excess density is called contrast, which is typically very small for biological objects containing light atoms only.

This is not true for neutron scattering where the variation of the contrast is a very powerful technique.

X-ray scattering

Scattering from dilute samples

Scattering from an ideal solution:

- no interaction between the particles
- only on species of particle
- particles are free to move

Dilute monodisperse solution fulfill these parameters!

X-ray scattering

Scattering from dilute samples

In the ideal dilute solutions the scattering intensity from the entire sample will be isotropic and proportional to the scattering from the single particle averaged over all orientations Ω .

Averaging of the term $\left\langle e^{\vec{s} \cdot \vec{r}_{jk}} \right\rangle_{\Omega} = \frac{\sin sr_{jk}}{sr_{jk}}$ lead to:

$$I(s) = \sum_j \sum_k \rho(r_j) \rho(r_k) \frac{\sin sr_{jk}}{sr_{jk}}$$

X-ray scattering

Scattering from dilute samples

The vectors s and r are now reduced to their absolute value, which lead to a significant loss of information in the SAXS pattern compared to e.g. crystallographic data

$$I(s) = \int_0^{\infty} p(r) \frac{\sin sr}{sr} dr$$

$s = |\vec{s}| = \frac{2\pi}{\lambda} \sin 2\Theta$


Debye 1915

The $p(r)$ function is called the pair distance distribution function.

The $I(s)$ is the Fourier-Transform of the $p(r)$ function.

X-ray scattering


The scattering function $I(s)$



X-Ray scattering

X-ray scattering

Pair distance distribution function $p(r)$




$$I(s) = 4\pi \int_{r=0}^{\infty} p(r) \frac{\sin(sr)}{sr} dr$$

$$p(r) = \frac{1}{2\pi^2} \int_{s=0}^{\infty} s^2 I(s) \frac{\sin(sr)}{sr} ds$$

Distance distribution function
"how often which distance" appears in the particle

X-ray scattering

Analysis of the scattering intensity



$$I(s) = \int_0^{\infty} p(r) \frac{\sin sr}{sr} dr$$

X-ray scattering

Radius of gyration – Guinier plot

Expanding the sin-function in the scattering intensity $I(s)$ in a series

$$\left(\frac{\sin sr}{sr} = 1 - \frac{s^2 r^2}{6} + \frac{s^4 r^4}{120} - \dots \right)$$

, the Debye-Function can be expressed as

$$I(s) \cong \int_0^{\infty} p(r) \left(1 - \frac{s^2 r^2}{6} + \dots \right) dr = \int_0^{\infty} p(r) dr - \frac{s^2}{6} \int_0^{\infty} p(r) r^2 dr + \dots$$

X-ray scattering

Radius of gyration – Guinier plot

The second moment of the distance distribution function $p(r)$ is given by:

$$\int_0^{\infty} p(r) r^2 dr$$

and the radius of gyration R_g can be introduced in analogy with classical mechanics:

$$R_g^2 = \frac{1}{2} \frac{\int_0^{\infty} p(r) r^2 dr}{\int_0^{\infty} p(r) dr}$$

X-ray scattering

Radius of gyration – Guinier plot

Putting the things together one gets:


$$I(s) \cong \int_0^{\infty} p(r) \left(1 - \frac{Rg^2}{6} s^2 + \dots \right) dr$$

And finally the Guinier approximation:

$$I(s) \cong I_0 e^{-\frac{Rg^2}{3} s^2}$$

Guinier and Fournet 1955

X-ray scattering Radius of gyration – Guinier plot



Guinier-Approximation:

- every scattering curve is for $sR_g < 1.2$ a Gauss function
- In the so called Guinier-Plot ($\ln I$ versus s^2) the R_g can be calculated from the slope of the plot
- The R_g is directly related to the form and the mass distribution of the particle

$$\ln I(s) = -\frac{R_g^2}{3} s^2 + \ln I_0$$

X-ray scattering




Radius of gyration – Guinier plot

The R_g is historically the first structural parameter determined from the SAXS data (Guinier and Fournet 1955), yielding direct information about the particle size and shape.


The value of I_0 is proportional to the squared number of electrons in the particle and to the particle concentration, so that absolute measurements allow one to determine the molecular mass of the solute.

X-Ray scattering


X-ray scattering Guinier approximation




Deviations from the ideal protein solution are easily visible in the Guinier plot!



Ideal solution of particles



Repulsive particle interactions




Attractive particle interactions

X-ray scattering Kratky plot

Kratky plots ($I(s) \cdot s^2$ versus s) can be used to identify disordered states and distinguish them from globular particles.

The scattering intensity of a globular protein has a Gaussian-like shape at small s and decays approximately as $1/s^4$ at high s yielding a bell-shaped Kratky plot with a well defined maximum.

In the case of an unfolded protein, the Kratky plot also presents a plateau over a specific range of s , which is followed by a monotonic increase.






X-Ray scattering

X-ray scattering Kratky plot

The Kratky plot is typically used to analyze the conformation of proteins, but can be used to analyze the random walk model of polymers.

A Kratky plot can be made by plotting:

$$I(s) \cdot s^2 \text{ versus } s.$$



Folded globular protein

Completely unfolded protein

Partially folded protein

X-ray scattering

Porod invariant

Another important overall value computed from the SAXS data is the so-called Porod invariant Q :

$$\int_0^{\infty} s^2 I(s) ds = Q$$

Porod 1951

Contrary to R_g , the Porod invariant depends only on the particle volume and not on its form.

For the Porod analysis, the behavior of the scattering intensity at higher scattering vectors plays a significant role.

This higher angle part of the $I(s)$ corresponds to the small interparticle distances. Therefore in the high s regime the scattering from the internal structure and in particular near the particle surface S is dominating the signal.

X-Ray scattering

X-ray scattering Porod volume

For compact particles, the asymptotic behavior at large s is described as

$$I(s) \approx \frac{K}{s^4}$$

introducing K as a constant depended on surface of the scattering particle.

This specific surface value can be estimated directly from the asymptotic behavior of the plot

$s^4 I(s)$ versus s

(Porod plot) at high s .

X-Ray scattering

X-ray scattering Porod plot

Porod approximation:


$$I(s) = K \cdot s^{-4}$$

$$\frac{K}{Q} \propto \frac{S}{V}$$

$$Q = \int_{s=0}^{\infty} s^2 I(s) ds$$

Q = Porod invariant

$$\frac{S}{V} \text{ shape to volume ratio}$$




Rule of thumb: Porod volume is approximately two times molecular weight

Here BSA: MW 66 kDa; Porod volume 115

Hardware and instrumentation part

Biological Small Angle X-ray Scattering:

1. Biological SAXS is a solution scattering technique
 - a. Change of environmental parameters such as temperature, pH, salt is possible
 - b. No protein crystals necessary
2. “Single shot” technique
 - a. Only 1-dim data recorded



Low resolution method
No direct recalculation of 3-dim model

X-Ray scattering

SAXS Synchrotron based

At all running synchrotron source are SAXS beamlines available!


At all new constructed synchrotrons SAXS beamlines are planned!

All SAXS beamlines are highly oversubscribed and “working horses” on their facilities!


The scientific applications ranging from soft condensed matter, nano-science and fiber diffraction on ordered biological systems up to structural biology in solution.




HASYLab, Hamburg




ESRF, Grenoble



APS, Chicago



SPring 8, Hyogo



LÜbeck 2009





EMBO Course 2012

X-Ray scattering

SAXS Lab sources




Bruker NanoStar



Bruker NanoStar
Anton Paar SAXSess
Rigaku BioSAXS 1000



Rigaku BioSAXS 1000




X-Ray scattering

Radiation from Synchrotron Storage Rings


Production of X-rays

Bending magnet




Necessary part of a synchrotron. Many of these dipole magnets form the synchrotron storage ring. The electrons (positrons) are deflected and accelerated in the magnetic field. This acceleration generates the synchrotron light. The light is emitted tangential to the electron beam.

Wiggler




Stack of magnetic dipoles. Put into the straight sections of the storage ring. Wigglers produce more light than bending magnets in a smaller source size.

Undulator




Most powerful insertion device! A stack of magnetic dipoles generate a high flux of photons in a very small source size. The specific arrangement of the dipoles ($d=n\lambda$) produces a discrete spectrum with coherent properties.



X-ray sources Brilliance




Increase of the beam
brilliance from X-ray tubes
to third generation
synchrotron sources to free
electron lasers (forth
generation)



X-Ray scattering

Radiation from Synchrotron Storage Rings

Dipole bending magnet (APS)



X-Ray scattering

X-ray optics Beam defining

Beam defining by slit pairs:




FIG. 2. Slit camera with three slits S_1 , S_2 and S_3 , each consisting of a pair of edges running perpendicular to the plane of paper. The dimensions in the vertical direction are greatly enlarged.

From Glatter & Kratky 1982

Small Angle Scattering


Book online available!

<http://physchem.kfunigraz.ac.at/sm/Index.html>

Or google "Otto Glatter"


Problem:

Every edge cutting the beam produces parasitic scattering! Refractive streaks at low angles impedes SAXS experiments!



Solution:

Successive slit systems. First slits cutting the beam and second system cuts out the undesired parasitic scattering



X-Ray scattering

X-ray Optics Monochromator


Monochromatization of the X-ray is achieved by using single crystals in Bragg diffraction geometry

$$n\lambda = 2d \sin(\Theta)$$

For instance a Si 111 crystal with d-spacing (=distance between the crystal lattice) of $d=3.14 \text{ \AA}$ deflects the beam for a wavelength of $\lambda=1 \text{ \AA}$ to $\Theta \sim 10^\circ$ (9.16°).

Problem:

The Bragg condition is also fulfilled for multiple orders of the wavelength n . These higher order have to be filtered by the X-ray mirror





X-Ray scattering

X-ray optics


Focusing

The glancing angle for X-rays on surfaces are very small. Only under grazing conditions X-ray mirrors can be used for focusing.

Bending of the highly polished mirror surface permits focusing on the parabolic mirror profile




Reflectivity for Rh@8keV



In order to use the full beam size X-ray mirrors on synchrotrons are typically in the range of 30cm to 1m. The grazing angle can be increased by using high Z-elements (e.g. Rhodium, Platin, Nickel) as reflecting surface. These mirrors are acting as well as high energy filter; important for monochromatization

Schematic X33 SAXS setup



The PETRA-III BioSAXS beamline P12




Detector tube and automated sample changer

X-Ray scattering



X-ray detectors

Single photon counting detectors



The PILATUS pixel detectors are at the moment the best detectors for solution SAXS.

Data acquisition integrated into the beamline control software (BMS). Advanced data storage concept for fast operation of a combined SAXS/WAXS detector setup.



X-Ray scattering

X-ray detectors Single photon counting detectors


Bump bonding over small indium metal balls.

Advanced microtechnology for production needed.

Every pixel is readout separately.

Single photon counting device.


Expensive technology!



X-Ray scattering

X-ray detectors

Low background noise recording



We are looking to protein signals of less than 0.5% above the back- ground level!

Summary and Conclusion

X-rays are scattered at the electrons of the atomic shell.

The heavier an atom the better it scatters! Proteins consist mostly of hydrogen, carbon and oxygen

For solution scattering the particles have to be in an ideal solution.

Due to the random orientation of the particles in solution, no atomic resolution is possible!

Summary and Conclusion

The Guinier approximation is valid for small s values, i.e. large interparticle distances. The radius of gyration is related to the shape of the particle

From the shape of the scattering function interparticle interactions can be estimated.


The Porod and Kratky analysis is valid at the high s regime and the scattering occurs from the internal structure. The Kratky plot gives information of the folding state of a protein.

From the Porod analysis the specific volume of a protein can be derived.


X-Ray scattering

X-ray scattering vs crystallography

Protein crystallography: Molecules (proteins) are fixed in a discrete crystal lattice




$$I(\vec{s}) = \int_V \rho^*(\vec{r})^2 e^{-i\vec{s}\vec{r}} dV$$



3-dim information about:
Length of the vector **AB**
Orientation of vector **AB**

Fourier
Transformation
→


3-dim protein structure
with atomic resolution



X-Ray scattering

X-ray scattering vs crystallography

Small angle scattering investigates protein solutions!




$$I(\vec{s}) = \int_V \rho^*(\vec{r})^2 e^{-i\vec{s}\vec{r}} dV$$

$$\langle e^{-i\vec{s}\vec{r}} \rangle_\Omega = \frac{\sin(sr)}{sr}$$


$$I(s) = 4\pi \int_0^\infty \rho^*(r) r^2 \frac{\sin(sr)}{sr} dr$$

$$I(s) = 4\pi \int_0^\infty p(r) \frac{\sin(sr)}{sr} dr$$

$$s = \frac{1}{d} = \frac{2\pi}{\lambda} \sin(2\Theta)$$




Radial
average



1-dim data

X-Ray scattering

Protein crystallography



Protein SAXS solution scattering



X-Ray scattering

SAXS is a cool technique!

You can use a lot of expensive and funny machines!

Let's meet at the beamline(s)!

